

Accelerating Research with AWS IoT

Marcilio Mendonca¹, Bruno Vitali¹, Arjun Shakhder¹, Brent Maranzano², Giuseppe Cogoni², Bo Du², Bryan Stone¹, and Don Bennett¹

¹AWS

²Pfizer, Inc.

The speed of designing synthetic routes and defining formulations for new drug substances can be limited by the rate at which the data becomes available for decision making. Much of the data is generated by offline laboratory experiments, which requires process sampling, subdividing, transporting and consequently time-consuming sample preparation prior to analysis, in addition to the time to analyze, interpret and report the results. An obvious path to accelerate sample measurements is to increase the number of analyses that can be performed in parallel by either increasing the staff or through automation. Alternatively, the entire process can be shortened by performing the analysis in-situ. This approach, referred to as Process Analytical Technology (PAT), has been used in the pharmaceutical industry for decades, but has been limited to monitoring a few key parameters, most commonly using spectroscopic tools. The advent of miniature, cheap electronics coupled with the Internet of Things (IoT) has opened the possibility of monitoring a wide range of parameters in real-time.

Motivations

The one obstacle that has inhibited the widespread adoption of IoT in the pharmaceutical industry is the lack of a secure, scalable and cost-effective platform that can be used to store, analyze and visualize the data.

In this post, it is described the development of a combination of on-premise and cloud-based IoT platform that can be used to monitor and analyze data from a wide range of sensors. This was achieved through the Pfizer's SmartLab initiative and in collaboration with AWS to integrate its IoT services.

Architecture

The key to a successful IoT platform is the ability to scale to a large number of devices that can be readily implemented in various processes with minimal experience or training by the end user e.g.: "plug-n-play". This minimal training interface often requires developing custom in-house applications that make calibration, measurement and analysis simple and robust. Moreover, the IoT platform must handle the large amount of data generated by the sensors securely and cost-effectively. One way to achieve this goal is through the use of AWS IoT services and solutions, which allow to connect, collect, store and analyze IoT data for industrial workloads.

At Pfizer, for this PoC, the PAT sensors communicate to the AWS IoT Core broker through a Mosquitto MQTT broker/bridge that communicates from the Pfizer lab network to a Rapid Virtual Private Cloud (VPC) sub-network. A web dashboard was also developed to be able to visualize and control each instrument that connects to a public AWS VPC and through a transit gateway, to the Rapid VPC for the data exchange. A schematic, summarizing this hybrid on-premise/cloud architecture is shown on Figure 1.




Figure 1: Pfizer/AWS IoT architecture.

Interfacing the laboratory instruments

For this PoC, three commonly used PAT sensors were utilized in order to cover different data types that can be exchanged, these include:

1. A level sensor, for process safety applications (generating Boolean data type);
2. A pH sensor (generating scalar data type);
3. An UV spectrometer (generating array data type).

Each of instrument is connected to a hardware wrapper (Raspberry Pi™ or ESP32 microcontroller) that aims to convert and exchange the control commands and data from its native protocol into a standard-based messaging protocol, MQTT, following the payload structure according to the Eclipse Foundation [Sparkplug™ B](#) specification. An example of data payload (DDATA) coming from one of the instrument is showed in Figure 2.

```
{
  "timestamp": "1701199950000",
  "seq": 14,
  "metrics": [
    {
      "name": "status",
      "timestamp": "1701199950000",
      "dataType": "UInt32",
      "value": 0
    },
    {
      "name": "level_high",
      "timestamp": "1701199950000",
      "dataType": "UInt32",
      "value": 0
    }
  ]
}
```

Figure 2: Example of MQTT payload following Sparkplug™B.

The local communication within the lab network is then managed through an Eclipse Mosquitto™ MQTT broker, configured to bridge and share the communication with the AWS cloud managed MQTT broker, [IoT Core](#) service, residing within a Rapid VPC, as a private sub-network.

Data storage visualization and instrument control

The data generated by each instrument is locally stored within the Pfizer lab network into an Influx database as well as within the AWS cloud, into a Dynamo DB, this last is then utilized as the data source for the web dashboard, used as instrument control platform and data visualization. The web dashboard interfaces to the Rapid VPC, through a transit gateway, that communicated to the public AWS VPC where it uses the Cognito service of AWS for user authentication.




Figure 3: Data visualization web app dashboard.

Once connected, the user has the option to visualize the data (Figure 3) and/or calibrate the instrument(s) (Figure 4). The instrument calibration was part of the PoC to test the bidirectionality of the communication and the ability to remotely control the individual instruments.

Conclusions

In this blog post, we have described how Pfizer through its SmartLab initiative, alongside Amazon AWS was able to implement and deploy a secure, scalable and cost-effective platform that can be used to store, analyze and visualize the data coming from PAT instruments.

With the discussed PoC it will be possible to modularize the deployment to multiple lab areas, also across multiple sites and be able to access to the data and control each individual instrument virtually from anywhere within the company.




Figure 4: Instrument calibration panel withing webapp.